

For Research Use Only

# ANGPTL4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67577-3-PBS



## Basic Information

Catalog Number:

67577-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG30070

GenBank Accession Number:

BC023647

GeneID (NCBI):

51129

UNIPROT ID:

Q9BY76

Full Name:

angiopoietin-like 4

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

2D9B9

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sandwich ELISA standard curve of MP50037-1, ANGPTL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67577-2-PBS. Detection antibody: 67577-3-PBS. Standard: Ag30070. Range: 0.313-20 ng/mL