

For Research Use Only

# SSBP1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67582-4-PBS



## Basic Information

Catalog Number:

67582-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29976

GenBank Accession Number:

BC000895

GeneID (NCBI):

6742

UNIPROT ID:

Q04837

Full Name:

single-stranded DNA binding protein

1

Calculated MW:

148 aa, 17 kDa

Purification Method:

Protein G purification

CloneNo.:

1A12H1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50750-2, SSBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67582-4-PBS. Detection antibody: 67582-5-PBS. Standard:Ag29976. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50750-3, SSBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67582-4-PBS. Detection antibody: 67582-6-PBS. Standard:Ag29976. Range: 0.098-3.125 ng/mL