

For Research Use Only

Bovine IgG Monoclonal antibody

Catalog Number: 67615-1-Ig



Basic Information

Catalog Number:

67615-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Observed MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

1C5C5

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

bovine

Positive Controls:

WB : FBS, BGG

Background Information

Bovine IgG is named as BGG, bovine IgG heavy chain. This antibody can recognize heavy chain of bovine IgG.

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



62.5 nL fetal bovine serum was subjected to SDS-PAGE followed by western blot with 67615-1-Ig (Bovine IgG antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.