

For Research Use Only

CKAP5 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67631-2-PBS

Basic Information

Catalog Number:

67631-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28770

GenBank Accession Number:

BC120869

GeneID (NCBI):

9793

UNIPROT ID:

Q14008

Full Name:

cytoskeleton associated protein 5

Purification Method:

Protein A purification

CloneNo.:

1E1D5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50224-1, CKAP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67631-2-PBS. Detection antibody: 67631-1-PBS. Standard:Ag28770. Range: 0.391-100 ng/mL.