

For Research Use Only

# PPM1B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67647-3-PBS



## Basic Information

Catalog Number: 67647-3-PBS	GenBank Accession Number: BC012002	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 5495	CloneNo.: 2B8E1
Source: Mouse	UNIPROT ID: O75688	
Isotype: IgG1	Full Name: protein phosphatase 1B (formerly 2C), magnesium-dependent, beta isoform	
Immunogen Catalog Number: AG3871	Calculated MW: 479 aa, 53 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

Selected Validation Data



Sample test of MP51322-2, PPM1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67647-1-PBS. Detection antibody: 67647-3-PBS.



Cytometric bead array standard curve of MP51322-2, PPM1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67647-1-PBS. Detection antibody: 67647-3-PBS. Standard:Ag3871. Range: 0.195-200 ng/mL.