

For Research Use Only

PPM1B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67647-4-PBS



Basic Information

Catalog Number: 67647-4-PBS	GenBank Accession Number: BC012002	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 5495	CloneNo.: 1G2B5
Source: Mouse	UNIPROT ID: O75688	
Isotype: IgG1	Full Name: protein phosphatase 1B (formerly 2C), magnesium-dependent, beta isoform	
Immunogen Catalog Number: AG3871	Calculated MW: 479 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51322-3, PPM1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67647-1-PBS. Detection antibody: 67647-4-PBS.



Cytometric bead array standard curve of MP51322-3, PPM1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67647-1-PBS. Detection antibody: 67647-4-PBS. Standard:Ag3871. Range: 0.195-200 ng/mL