

For Research Use Only

# PPM1B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67647-5-PBS



## Basic Information

Catalog Number:

67647-5-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3871

GenBank Accession Number:

BC012002

GeneID (NCBI):

5495

UNIPROT ID:

O75688

Full Name:

protein phosphatase 1B (formerly 2C),  
magnesium-dependent, beta isoform

Calculated MW:

479 aa, 53 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A9A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sample test of MP51322-4, PPM1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67647-1-PBS. Detection antibody: 67647-5-PBS.



Cytometric bead array standard curve of MP51322-4, PPM1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67647-1-PBS. Detection antibody: 67647-5-PBS. Standard:Ag3871. Range: 0.195-200 ng/mL