

For Research Use Only

EPCR/CD201 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67658-4-PBS



Basic Information

Catalog Number:

67658-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG25628

GenBank Accession Number:

BC014451

GeneID (NCBI):

10544

UNIPROT ID:

Q9UNN8

Full Name:

protein C receptor, endothelial (EPCR)

Calculated MW:

238 aa, 27 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1C3A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50576-2, PROCR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67658-4-PBS. Detection antibody: 67658-5-PBS. Standard:Ag25628. Range: 1.563-200 ng/mL.