

For Research Use Only

RAB18 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67659-2-PBS



Basic Information

Catalog Number:

67659-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16669

GenBank Accession Number:

BC015014

GeneID (NCBI):

22931

UNIPROT ID:

Q9NP72

Full Name:

RAB18, member RAS oncogene family

Calculated MW:

206 aa, 23 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B9D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51469-1, RAB18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67659-1-PBS. Detection antibody: 67659-2-PBS.



Cytometric bead array standard curve of MP51469-1, RAB18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67659-1-PBS. Detection antibody: 67659-2-PBS. Standard: Ag16669. Range: 3.125-400 ng/mL.