

For Research Use Only

SNX5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67665-2-PBS



Basic Information

Catalog Number:

67665-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12330

GenBank Accession Number:

BC093623

GeneID (NCBI):

27131

UNIPROT ID:

Q9Y5X3

Full Name:

sorting nexin 5

Calculated MW:

404 aa, 47 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A2A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51693-1, SNX5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67665-2-PBS. Detection antibody: 67665-3-PBS.



Cytometric bead array standard curve of MP51693-1, SNX5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67665-2-PBS. Detection antibody: 67665-3-PBS. Standard:Ag12330. Range: 0.391-50 ng/mL.