

For Research Use Only

ALPK1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67689-4-PBS

Basic Information

Catalog Number:

67689-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5320

GenBank Accession Number:

BC060780

GeneID (NCBI):

80216

UNIPROT ID:

Q96QP1

Full Name:

alpha-kinase 1

Calculated MW:

139 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G3B4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50186-1, ALPK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67689-3-PBS. Detection antibody: 67689-4-PBS. Standard:Ag5320. Range: 1.563-100 ng/mL.