

For Research Use Only

# LAMB1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67705-3-PBS



## Basic Information

Catalog Number:

67705-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19182

GenBank Accession Number:

BC113455

GeneID (NCBI):

3912

UNIPROT ID:

P07942

Full Name:

laminin, beta 1

Calculated MW:

1786 aa, 198 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C1B4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

Cytometric bead array standard curve of MP51088-1, LAMB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67705-2-PBS. Detection antibody: 67705-3-PBS. Standard:null. Range: null.



Cytometric bead array standard curve of MP51088-1, LAMB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67705-2-PBS. Detection antibody: 67705-3-PBS. Standard:Ag19182. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51088-2, LAMB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67705-4-PBS. Detection antibody: 67705-3-PBS. Standard:Ag19182. Range: 0.195-50 ng/mL