

For Research Use Only

CEP290 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67716-3-PBS

Basic Information

Catalog Number: 67716-3-PBS	GenBank Accession Number: BC008641	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 80184	CloneNo.: 2H4G5
Source: Mouse	UNIPROT ID: O15078	
Isotype: IgG1	Full Name: centrosomal protein 290kDa	
Immunogen Catalog Number: AG17457	Calculated MW: 2479 aa, 290 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50156-1, CEP290 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67716-2-PBS. Detection antibody: 67716-3-PBS. Standard:Ag17457. Range: 1.563-100 ng/mL.