

# GIMAP2 Monoclonal antibody

Catalog Number: 67720-1-Ig

## Basic Information

Catalog Number:

67720-1-Ig

Size:

1500 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4368

GenBank Accession Number:

BC032345

GeneID (NCBI):

26157

UNIPROT ID:

Q9UG22

Full Name:

GTPase, IMAP family member 2

Calculated MW:

337 aa, 38 kDa

Observed MW:

33-38 kDa

Purification Method:

Protein G purification

CloneNo.:

1F2A2

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : MOLT-4 cells, Jurkat cells

## Background Information

### Storage

Storage:


Store at -20°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

## Selected Validation Data



MOLT-4 cells were subjected to SDS PAGE followed by western blot with 67720-1-Ig (GIMAP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.