

For Research Use Only

# CTPS2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67734-2-PBS



## Basic Information

Catalog Number:

67734-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30258

GenBank Accession Number:

BC034986

GeneID (NCBI):

56474

UNIPROT ID:

Q9NRF8

Full Name:

CTP synthase II

Calculated MW:

586 aa, 66 kDa

Purification Method:

Protein G purification

CloneNo.:

1D10F9

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sandwich ELISA standard curve of MP50097-1, CTPS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67734-1-PBS. Detection antibody: 67734-2-PBS. Standard:Ag30258. Range: 0.195-50 ng/mL