

For Research Use Only



# Alpha smooth muscle actin specific Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67735-2-PBS

## Basic Information

Catalog Number:

67735-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM\_001613

GeneID (NCBI):

59

UNIPROT ID:

P62736

Full Name:

actin, alpha 2, smooth muscle, aorta

Calculated MW:

42 kDa

Purification Method:

Protein G purification

CloneNo.:

1D7H12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50296-2, ACTA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68895-2-PBS. Detection antibody: 67735-2-PBS. Standard:Ag6088. Range: 3.125-50 ng/mL