

For Research Use Only

CCDC45 Monoclonal antibody, PBS Only



Catalog Number: 67755-1-PBS

Basic Information

Catalog Number:

67755-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27836

GenBank Accession Number:

BC009518

GeneID (NCBI):

90799

UNIPROT ID:

Q96GE4

Full Name:

coiled-coil domain containing 45

Calculated MW:

821 aa, 95 kDa

Observed MW:

95-100 kDa

Purification Method:

Protein G purification

CloneNo.:

1A12A7

Applications

Tested Applications:

WB, Indirect ELISA, IF

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67755-1-Ig (CCDC45 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67755-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed U2OS cells using CCDC45 antibody (67755-1-Ig, Clone: 1A12A7) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67755-1-PBS in a different storage buffer formulation.