

For Research Use Only

GPX4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67763-2-PBS



Basic Information

Catalog Number:

67763-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG30650

GenBank Accession Number:

BC021567

GeneID (NCBI):

2879

UNIPROT ID:

P36969

Full Name:

glutathione peroxidase 4
(phospholipid hydroperoxidase)

Purification Method:

Protein A purification

CloneNo.:

1A1A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50524-1, GPX4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67763-2-PBS. Detection antibody: 67763-3-PBS. Standard:Ag30650. Range: 0.098-12.5 ng/mL