

For Research Use Only

ALCAM Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67768-2-PBS



Basic Information

Catalog Number: 67768-2-PBS	GenBank Accession Number: BC137097	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 214	CloneNo.: 1C5H8
Isotype: IgG1	UNIPROT ID: Q13740	
Immunogen Catalog Number: AG16892	Full Name: activated leukocyte cell adhesion molecule	
	Calculated MW: 583 aa, 65 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50574-1, ALCAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67768-2-PBS. Detection antibody: 67768-3-PBS. Standard:Ag16892. Range: 3.125-100 ng/mL.