

For Research Use Only

PSMD2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67770-2-PBS



Basic Information

Catalog Number: 67770-2-PBS	GenBank Accession Number: BC002368	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 5708	CloneNo.: 1D12A9
Isotype: IgG2a	UNIPROT ID: Q13200	
Immunogen Catalog Number: AG6543	Full Name: proteasome (prosome, macropain) 26S subunit, non-ATPase, 2	
	Calculated MW: 100 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51750-1, PSMD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67770-2-PBS. Detection antibody: 67770-3-PBS.



Cytometric bead array standard curve of MP51750-1, PSMD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67770-2-PBS. Detection antibody: 67770-3-PBS. Standard:Ag6543. Range: 0.781-100 ng/mL.