

For Research Use Only

PSMD2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67770-3-PBS



Basic Information

Catalog Number:

67770-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6543

GenBank Accession Number:

BC002368

GeneID (NCBI):

5708

UNIPROT ID:

Q13200

Full Name:

proteasome (prosome, macropain)
26S subunit, non-ATPase, 2

Calculated MW:

100 kDa

Purification Method:

Protein G purification

CloneNo.:

3A1H10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51750-1, PSMD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67770-2-PBS. Detection antibody: 67770-3-PBS.



Cytometric bead array standard curve of MP51750-1, PSMD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67770-2-PBS. Detection antibody: 67770-3-PBS. Standard:Ag6543. Range: 0.781-100 ng/mL.