

For Research Use Only

# GALE Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67772-2-PBS



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
67772-2-PBS	BC050685	Protein A purification
Concentration:	GeneID (NCBI):	CloneNo.:
1 mg/ml	2582	3F1B10
Source:	UNIPROT ID:	
Mouse	Q14376	
Isotype:	Full Name:	
IgG2a	UDP-galactose-4-epimerase	
Immunogen Catalog Number:	Calculated MW:	
AG6571	38 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51281-1, GALE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67772-2-PBS. Detection antibody: 67772-3-PBS. Standard:Ag6571. Range: 6.25-100 ng/mL