

For Research Use Only

RAB11B Monoclonal antibody

Catalog Number: 67780-1-Ig



Basic Information

Catalog Number:

67780-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25429

GenBank Accession Number:

BC110081

GeneID (NCBI):

9230

UNIPROT ID:

Q15907

Full Name:

RAB11B, member RAS oncogene family

Calculated MW:

24 kDa

Observed MW:

24 kDa

Purification Method:

Protein G purification

CloneNo.:

1G7D6

Recommended Dilutions:

WB 1:5000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : LNCaP cells, NIH/3T3 cells, HSC-T6 cells, K-562 cells, Jurkat cells, HepG2 cells, HEK-293 cells, HeLa cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67780-1-Ig (RAB11B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control.