

For Research Use Only

RAB11B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67780-2-PBS



Basic Information

Catalog Number:

67780-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25429

GenBank Accession Number:

BC110081

GeneID (NCBI):

9230

UNIPROT ID:

Q15907

Full Name:

RAB11B, member RAS oncogene family

Calculated MW:

24 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D2B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51735-1, RAB11B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67780-1-PBS. Detection antibody: 67780-2-PBS.



Cytometric bead array standard curve of MP51735-1, RAB11B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67780-1-PBS. Detection antibody: 67780-2-PBS. Standard: Ag25429. Range: 1.563-200 ng/mL