

For Research Use Only

CS Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67784-3-PBS



Basic Information

Catalog Number:

67784-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG9255

GenBank Accession Number:

BC010106

GeneID (NCBI):

1431

UNIPROT ID:

O75390

Full Name:

citrate synthase

Calculated MW:

466 aa, 52 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2D7B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50606-1, CS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67784-2-PBS. Detection antibody: 67784-3-PBS. Standard:Ag9255. Range: 0.098-50 ng/mL



Cytometric bead array standard curve of MP50606-3, CS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67784-3-PBS. Detection antibody: 67784-5-PBS. Standard:Ag9255. Range: 0.098-100 ng/mL