

FGFR4 Monoclonal antibody

Catalog Number: 67800-1-Ig **1 Publications**

Basic Information

Catalog Number: 67800-1-Ig	GenBank Accession Number: BC011847	Purification Method: Protein G purification
Size: 1000 µg/ml	GeneID (NCBI): 2264	CloneNo.: 1A2B5
Source: Mouse	UNIPROT ID: P22455	Recommended Dilutions: WB 1:1000-1:6000
Isotype: IgG1	Full Name: fibroblast growth factor receptor 4	
Immunogen Catalog Number: AG17041	Calculated MW: 88 kDa	
	Observed MW: 100-110 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Raji cells, MOLT-4 cells, Jurkat cells
Cited Applications: WB	
Species Specificity: Human	
Cited Species: human	

Background Information


Notable Publications

Author	Pubmed ID	Journal	Application
Nianli Liu	37817227	J Exp Clin Cancer Res	WB

Storage

Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67800-1-Ig (FGFR4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.