

For Research Use Only

FGFR4 Monoclonal antibody

Catalog Number: 67800-1-Ig **1 Publications**



Basic Information

Catalog Number: 67800-1-Ig	GenBank Accession Number: BC011847	Purification Method: Protein G purification
Size: 1000 µg/ml	GeneID (NCBI): 2264	CloneNo.: 1A2B5
Source: Mouse	UNIPROT ID: P22455	Recommended Dilutions: WB 1:1000-1:6000
Isotype: IgG1	Full Name: fibroblast growth factor receptor 4	
Immunogen Catalog Number: AG17041	Calculated MW: 88 kDa	
	Observed MW: 100-110 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Raji cells, MOLT-4 cells, Jurkat cells
Cited Applications: WB	
Species Specificity: Human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Nianli Liu	37817227	J Exp Clin Cancer Res	WB


Storage

Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67800-1-Ig (FGFR4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.