

For Research Use Only

SLC25A42 Monoclonal antibody

Catalog Number: 67819-1-Ig



Basic Information

Catalog Number:

67819-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG26379

GenBank Accession Number:

BC045598

GeneID (NCBI):

284439

UNIPROT ID:

Q86VD7

Full Name:

solute carrier family 25, member 42

Calculated MW:

318 aa, 35 kDa

Observed MW:

33-35 kDa

Purification Method:

Protein A purification

CloneNo.:

2F7G10

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : K-562 cells, LNCaP cells, Jurkat cells, Raji cells

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 67819-1-Ig (SLC25A42 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.