

For Research Use Only

# PPARGC1B Monoclonal antibody, PBS Only



Catalog Number: 67821-1-PBS

## Basic Information

Catalog Number:

67821-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28861

GenBank Accession Number:

BC132971

GeneID (NCBI):

133522

UNIPROT ID:

Q86YN6

Full Name:

peroxisome proliferator-activated  
receptor gamma, coactivator 1 beta

Calculated MW:

1023 aa, 113 kDa

Observed MW:

150 kDa

Purification Method:

Protein G purification

CloneNo.:

2D2E8

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, pig, rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67821-1-Ig (PPARGC1B antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67821-1-PBS in a different storage buffer formulation.



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 67821-1-Ig (PPARGC1B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67821-1-PBS in a different storage buffer formulation.