

For Research Use Only

Caspase 6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67825-3-PBS



Basic Information

Catalog Number:

67825-3-PBS

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG20534

GenBank Accession Number:

BC000305

GeneID (NCBI):

839

UNIPROT ID:

P55212

Full Name:

caspase 6, apoptosis-related cysteine peptidase

Calculated MW:

33 kDa, 22 kDa

Purification Method:

Protein A purification

CloneNo.:

2D1C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51728-1, Caspase 6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67825-2-PBS. Detection antibody: 67825-3-PBS.



Cytometric bead array standard curve of MP51728-1, Caspase 6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67825-2-PBS. Detection antibody: 67825-3-PBS. Standard:Ag20534. Range: 6.25-200 ng/mL