

For Research Use Only

ID1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67827-3-PBS



Basic Information

Catalog Number:

67827-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13359

GenBank Accession Number:

BC000613

GeneID (NCBI):

3397

UNIPROT ID:

P41134

Full Name:

inhibitor of DNA binding 1, dominant negative helix-loop-helix protein

Calculated MW:

16 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H4D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51703-1, ID1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67827-3-PBS. Detection antibody: 67827-2-PBS.



Cytometric bead array standard curve of MP51703-1, ID1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67827-3-PBS. Detection antibody: 67827-2-PBS. Standard: Ag13359. Range: 0.391-50 ng/mL.