

For Research Use Only

INADL Monoclonal antibody, PBS Only



Catalog Number: 67833-1-PBS

Basic Information

Catalog Number: 67833-1-PBS	GenBank Accession Number: BC142661	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 10207	CloneNo.: 3C6A4
Source: Mouse	UNIPROT ID: Q8NI35	
Isotype: IgG1	Full Name: InaD-like (Drosophila)	
Immunogen Catalog Number: AG15953	Calculated MW: 1801 aa, 196 kDa	
	Observed MW: 196 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67833-1-Ig (INADL antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67833-1-PBS in a different storage buffer formulation.