

For Research Use Only

ELAVL4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67835-3-PBS



Basic Information

Catalog Number:

67835-3-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG3678

GenBank Accession Number:

BC036071

GeneID (NCBI):

1996

UNIPROT ID:

P26378

Full Name:

ELAV (embryonic lethal, abnormal vision, Drosophila)-like 4 (Hu antigen D)

Calculated MW:

380 aa, 42 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2B1B6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51680-1, ELAVL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67835-2-PBS. Detection antibody: 67835-3-PBS.



Cytometric bead array standard curve of MP51680-1, ELAVL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67835-2-PBS. Detection antibody: 67835-3-PBS. Standard:Ag3678. Range: 0.781-100 ng/mL