

For Research Use Only

CAPNS1 Monoclonal antibody

Catalog Number: 67838-1-Ig



Basic Information

Catalog Number: 67838-1-Ig	GenBank Accession Number: BC011903	Purification Method: Protein G purification
Size: 1000 µg/ml	GeneID (NCBI): 826	CloneNo.: 2B6G6
Source: Mouse	UNIPROT ID: P04632	Recommended Dilutions: WB 1:2000-1:10000 IF 1:250-1:1000
Isotype: IgG1	Full Name: calpain, small subunit 1	
Immunogen Catalog Number: AG21885	Calculated MW: 28 kDa	

Applications

Tested Applications:
IF/ICC, WB, ELISA

Species Specificity:
Human, rat

Positive Controls:

WB: A375 cells, A431 cells, HeLa cells, K-562 cells, THP-1 cells, HSC-T6 cells, PC-12 cells

IF: A431 cells,

Background Information

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A375 cells were subjected to SDS PAGE followed by western blot with 67838-1-Ig (CAPNS1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using CAPNS1 antibody (67838-1-Ig, Clone: 2B6G6) at dilution of 1:500 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).