

For Research Use Only

CAPNS1 Monoclonal antibody

Catalog Number: 67838-1-Ig



Basic Information

Catalog Number:

67838-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG21885

GenBank Accession Number:

BC011903

GeneID (NCBI):

826

UNIPROT ID:

P04632

Full Name:

calpain, small subunit 1

Calculated MW:

28 kDa

Purification Method:

Protein G purification

CloneNo.:

2B6G6

Recommended Dilutions:

WB 1:2000-1:10000

IF/ICC 1:250-1:1000

Applications

Tested Applications:

IF/ICC, WB, ELISA

Species Specificity:

Human, rat

Positive Controls:

WB : A375 cells, A431 cells, HeLa cells, K-562 cells, THP-1 cells, HSC-T6 cells, PC-12 cells

IF/ICC : A431 cells,

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A375 cells were subjected to SDS PAGE followed by western blot with 67838-1-Ig (CAPNS1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using CAPNS1 antibody (67838-1-Ig, Clone: 2B6G6) at dilution of 1:500 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).