

For Research Use Only

PGA5 Monoclonal antibody

Catalog Number: 67848-1-Ig



Basic Information

Catalog Number:

67848-1-Ig

Size:

2287 μ g/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11657

GenBank Accession Number:

BC029055

GeneID (NCBI):

5222

UNIPROT ID:

P0DJ09

Full Name:

pepsinogen 5, group I (pepsinogen A)

Calculated MW:

388 aa, 42 kDa

Purification Method:

Protein G purification

CloneNo.:

1C5E12

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

