

For Research Use Only

# PGA5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67848-3-PBS



## Basic Information

Catalog Number:

67848-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11657

GenBank Accession Number:

BC029055

GeneID (NCBI):

5222

UNIPROT ID:

P0DJ09

Full Name:

pepsinogen 5, group I (pepsinogen A)

Calculated MW:

388 aa, 42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C11H3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sample test of MP51334-1, PGA5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67848-2-PBS. Detection antibody: 67848-3-PBS.



Cytometric bead array standard curve of MP51334-1, PGA5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67848-2-PBS. Detection antibody: 67848-3-PBS. Standard:Ag11657. Range: 3.125-50 ng/mL.