

For Research Use Only

# CD42b Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67857-3-PBS



## Basic Information

Catalog Number:

67857-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31181

GenBank Accession Number:

BC027955

GeneID (NCBI):

2811

ENSEMBL Gene ID:

ENSG00000185245

UNIPROT ID:

P07359

Full Name:

glycoprotein Ib (platelet), alpha polypeptide

Calculated MW:

626 aa, 69 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F10D3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51744-1, GP1BA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67857-2-PBS. Detection antibody: 67857-3-PBS.



Cytometric bead array standard curve of MP51744-1, GP1BA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67857-2-PBS. Detection antibody: 67857-3-PBS. Standard:Ag31181. Range: 0.781-100 ng/mL.