

For Research Use Only

AES Monoclonal antibody, PBS Only

Catalog Number: 67869-1-PBS



Basic Information

Catalog Number:

67869-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31320

GenBank Accession Number:

BC113735

GeneID (NCBI):

166

UNIPROT ID:

Q08117

Full Name:

amino-terminal enhancer of split

Calculated MW:

264 aa, 29 kDa

Observed MW:

31 kDa

Purification Method:

Protein G purification

CloneNo.:

1F7E1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67869-1-Ig (AES antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67869-1-PBS in a different storage buffer formulation.