

For Research Use Only

OSBPL5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67913-3-PBS



Basic Information

Catalog Number:

67913-3-PBS

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG11554

GenBank Accession Number:

BC032646

GeneID (NCBI):

114879

UNIPROT ID:

Q9H0X9

Full Name:

oxysterol binding protein-like 5

Calculated MW:

811 aa, 91 kDa

Purification Method:

Protein A purification

CloneNo.:

1B11E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50959-1, OSBP.L5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67913-2-PBS. Detection antibody: 67913-3-PBS. Standard:Ag11554. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 67913-3-Ig against Human RUVBL2 was performed. The affinity constant is below 1 pM.