

For Research Use Only

# SLC36A3 Monoclonal antibody, PBS Only (Capture)

Catalog Number:67929-2-PBS



## Basic Information

<b>Catalog Number:</b> 67929-2-PBS	<b>GenBank Accession Number:</b> BC101092	<b>Purification Method:</b> Protein G purification
<b>Source:</b> Mouse	<b>GeneID (NCBI):</b> 285641	<b>CloneNo.:</b> 2C1D1
<b>Isotype:</b> IgG1	<b>UNIPROT ID:</b> Q495N2	
<b>Immunogen Catalog Number:</b> AG16556	<b>Full Name:</b> solute carrier family 36 (proton/amino acid symporter), member 3	
	<b>Calculated MW:</b> 511 aa, 56 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

