

For Research Use Only

# SLC36A3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67929-3-PBS



## Basic Information

Catalog Number:

67929-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16556

GenBank Accession Number:

BC101092

GeneID (NCBI):

285641

UNIPROT ID:

Q495N2

Full Name:

solute carrier family 36  
(proton/amino acid symporter),  
member 3

Calculated MW:

511 aa, 56 kDa

Purification Method:

Protein G purification

CloneNo.:

1H9D7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

