

For Research Use Only

SCAMP3 Monoclonal antibody, PBS Only



Catalog Number: 67932-1-PBS

Basic Information

Catalog Number:

67932-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26073

GenBank Accession Number:

BC000161

GeneID (NCBI):

10067

UNIPROT ID:

O14828

Full Name:

secretory carrier membrane protein 3

Calculated MW:

38 kDa

Observed MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

1C1E6

Applications

Tested Applications:

WB, IF, ELISA

Species Specificity:

Human, rat, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67932-1-Ig (SCAMP3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67932-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using SCAMP3 antibody (67932-1-Ig, Clone: 1C1E6) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67932-1-PBS in a different storage buffer formulation.