

For Research Use Only

# WDR1 Monoclonal antibody

Catalog Number: 67937-1-Ig

Featured Product

1 Publications



## Basic Information

Catalog Number:

67937-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30967

GenBank Accession Number:

BC030541

GeneID (NCBI):

9948

UNIPROT ID:

O75083

Full Name:

WD repeat domain 1

Calculated MW:

606 aa, 66 kDa

Observed MW:

65-66 kDa

Purification Method:

Protein G purification

CloneNo.:

2D11E1

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, rat

Cited Species:

human

Positive Controls:

WB: A549 cells, rat testis tissue, HEK-293 cells, U-87 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Jingyun Li	37366253	Acta Physiol (Oxf)	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67937-1-Ig (WDR1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



WB result of WDR1 antibody (67937-1-Ig; 1:10000; incubated at room temperature for 1.5 hours) with sh-Control and sh-WDR1 transfected HEK-293 cells.