

For Research Use Only

PUF-A, KIAA0020 Monoclonal antibody



Catalog Number: 67941-1-Ig

Basic Information

Catalog Number:

67941-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31324

GenBank Accession Number:

BC016137

GeneID (NCBI):

9933

UNIPROT ID:

Q15397

Full Name:

KIAA0020

Calculated MW:

648 aa, 74 kDa

Observed MW:

68-70 kDa

Purification Method:

Protein G purification

CloneNo.:

1G4B7

Recommended Dilutions:

WB 1:2000-1:20000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : HeLa cells, HSC-T6 cells, PC-12 cells, NIH/3T3 cells, 4T1 cells, HEK-293 cells, HepG2 cells, Jurkat cells, k-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67941-1-Ig (KIAA0020 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control.