

For Research Use Only

ACF1 Monoclonal antibody, PBS Only

Catalog Number: 67947-1-PBS



Basic Information

Catalog Number:

67947-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG30607

GenBank Accession Number:

NM_013448

GeneID (NCBI):

11177

UNIPROT ID:

Q9NRL2

Full Name:

bromodomain adjacent to zinc finger domain, 1A

Calculated MW:

179 kDa

Purification Method:

Protein A purification

CloneNo.:

1C4F8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67947-1-Ig (ACF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67947-1-PBS in a different storage buffer formulation.