

For Research Use Only

Caspase 7 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67956-2-PBS



Basic Information

Catalog Number:	67956-2-PBS	GenBank Accession Number:	BC015799	Purification Method:	Protein G purification
Source:	Mouse	GeneID (NCBI):	840	CloneNo.:	1G9C1
Isotype:	IgG1	UNIPROT ID:	P55210		
Immunogen Catalog Number:	AG27601	Full Name:	caspase 7, apoptosis-related cysteine peptidase		
		Calculated MW:	303 aa, 34 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51724-1, Caspase 7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67956-2-PBS. Detection antibody: 67956-3-PBS.



Cytometric bead array standard curve of MP51724-1, Caspase 7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67956-2-PBS. Detection antibody: 67956-3-PBS. Standard:Ag27601. Range: 1.563-200 ng/mL.