

For Research Use Only

Caspase 7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67956-3-PBS



Basic Information

Catalog Number:

67956-3-PBS

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG27601

GenBank Accession Number:

BC015799

GeneID (NCBI):

840

UNIPROT ID:

P55210

Full Name:

caspase 7, apoptosis-related cysteine peptidase

Calculated MW:

303 aa, 34 kDa

Purification Method:

Euglobulin precipitation

CloneNo.:

1C1E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51724-1, Caspase 7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67956-2-PBS. Detection antibody: 67956-3-PBS.



Cytometric bead array standard curve of MP51724-1, Caspase 7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67956-2-PBS. Detection antibody: 67956-3-PBS. Standard:Ag27601. Range: 1.563-200 ng/mL.