

For Research Use Only

SATB2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67958-1-PBS

Featured Product



Basic Information

Catalog Number:

67958-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17237

GenBank Accession Number:

BC098136

GeneID (NCBI):

23314

UNIPROT ID:

Q9UPW6

Full Name:

SATB homeobox 2

Calculated MW:

733 aa, 83 kDa

Observed MW:

85-100 kDa

Purification Method:

Protein G purification

CloneNo.:

2F8E2

Applications

Tested Applications:

WB, IHC, Cytometric bead array, Indirect ELISA

Species Specificity:

human, mouse

Background Information

SATB2, also named as KIAA1034, belongs to the CUT homeobox family. SATB2 binds to DNA at nuclear matrix- or scaffold-associated regions. SATB2 recognizes the sugar-phosphate structure of double-stranded DNA. SATB2 is a transcription factor controlling nuclear gene expression, by binding to matrix attachment regions (MARs) of DNA and inducing a local chromatin-loop remodeling. SATB2 acts as a docking site for several chromatin remodeling enzymes and also by recruiting corepressors (HDACs) or coactivators (HATs) directly to promoters and enhancers. It is required for the initiation of the upper-layer neurons (UL1) specific genetic program and for the inactivation of deep-layer neurons (DL) and UL2 specific genes, probably by modulating BCL11B expression. It is a repressor of Ctip2 and regulatory determinant of corticocortical connections in the developing cerebral cortex. SATB2 may play an important role in palate formation. SATB2 acts as a molecular node in a transcriptional network regulating skeletal development and osteoblast differentiation. SATB2 has two isoforms produced by alternative splicing with the MW of 70 kDa and 83 kDa. It can be detected as 85-105 kDa by sumo modification (PMID: 14701874, PMID: 35140581).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 67958-1-Ig (SATB2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67958-1-PBS in a different storage buffer formulation.




WB result of SATB2 antibody (67958-1-Ig; 1:5000; incubated at room temperature for 1.5 hours) with sh-Control and sh-SATB2 transfected HEK-293 cells. This data was developed using the same antibody clone with 67958-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 67958-1-Ig (SATB2 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67958-1-PBS in a different storage buffer formulation.




HepG2 cell lysates and HepG2 nuclear protein were subjected to SDS PAGE followed by western blot with 67958-1-Ig (SATB2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67958-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer slide using 67958-1-Ig (SATB2 antibody) at dilution of 1:1000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67958-1-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP50682-1, SATB2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67958-1-PBS. Detection antibody: 67958-2-PBS. Standard: Ag17237. Range: 0.781-100 ng/mL.



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 67958-1-Ig (SATB2 antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67958-1-PBS in a different storage buffer formulation.