

For Research Use Only

# NFIX Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67983-2-PBS



## Basic Information

<b>Catalog Number:</b> 67983-2-PBS	<b>GenBank Accession Number:</b> NM_001271043	<b>Purification Method:</b> Protein G Magarose purification
<b>Source:</b> Mouse	<b>GeneID (NCBI):</b> 4784	<b>CloneNo.:</b> 1B2F5
<b>Isotype:</b> IgG1	<b>UNIPROT ID:</b> Q14938	
<b>Immunogen Catalog Number:</b> AG31447	<b>Full Name:</b> nuclear factor I/X (CCAAT-binding transcription factor)	
	<b>Calculated MW:</b> 50 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51754-1, NFIX Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67983-2-PBS. Detection antibody: 67983-1-PBS.



Cytometric bead array standard curve of MP51754-1, NFIX Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67983-2-PBS. Detection antibody: 67983-1-PBS. Standard: Ag31447. Range: 0.313-40 ng/mL.